

ISSI LC-PZA2-E

Canon® Ethernet Lens
Controller

Application Program Interface

Programmers Reference Manual

Copyright © ISSI June 2021

Version 4.7.19

Notice

Copyright © 2021 Innovative Scientific Solutions Inc (ISSI). All rights reserved.

ISSI does not guarantee that the ISSI LC-PZA2-E API will function properly in every hardware/software environment. This software is inherently complex, and users are cautioned to verify the results of their work. ISSI has tested the software and reviewed the documentation. ISSI MAKES NO WARRANTY OR REPRESENTATION, EITHER EXPRESSED OR IMPLIED, WITH RESPECT TO THIS SOFTWARE OR DOCUMENTATION, THEIR QUALITY, PERFORMANCE, MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE. AS A RESULT, THIS SOFTWARE AND DOCUMENTATION ARE LICENSED “AS IS” AND YOU, THE LICENSEE ARE ASSUMING THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE. IN NO EVENT WILL ISSI BE LIABLE FOR DIRECT, INDIRECT, SPECIAL, OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE SOFTWARE OR DOCUMENTATION, even if advised of the possibility of such damages. In particular, API shall have no liability for any programs or data stored or used with ISSI software, including the costs of recovering such programs or data. Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

1. Redistributions of source code must retain the above copyright notice, this list of conditions and the above disclaimer.
2. Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the above disclaimer listed in this license in the documentation and/or other materials provided with the distribution.
3. Neither the name of the copyright holders nor the names of its contributors may be used to endorse or promote products derived from this API without specific prior written permission.

Contact Information: Innovative Scientific Solutions Inc.

7610 McEwen Road Dayton, OH 45459 Ph.: (937) 630-3012 Fax: (937) 630-3015

Website: www.innssi.com, Email : issi-sales@innssi.com

ISSI LC-PZA2-E API

API based on UDP Protocol. Port 1339/UDP.

1. Commands

Action: Get controller f/w version

Command: "ver<CR>"

Returns: "ISSI LCI-EF-232=4.x.y (s/n:19xxx)"

Action: Initialize lens, and get current motors values and ranges (aperture ranges would change when zoom changed)

Command: "ping<CR>"

Returns (in separate UDP packets):

zRange=minZ,maxZ

fRange=minF,maxF

aRange=minA,maxA

Current=ZZ,FF,AA

AF=X

IS=Y

ISactive=S"

where

- minZ– minimum zoom value,
maxZ – maximum zoom value
- minF– minimum focus value,
maxF – maximum focus value

- minA– minimum aperture value,
maxA – maximum aperture value
- ZZ– current zoom value, FF– current focus value, AA - current aperture value
- X = **1** if the focus switch is in *AF* position on lens and **0** if in *MF*.
- Y = **1** if Image Stabilization(IS) switch on the lens is enabled, and **0** if disabled or the lens has no IS
- S = **1** if IS is activated and **0** if it deactivated
- “*nolensfound*” return in case no lens found

Action: Move Focus motor on X units

Command: “*moveFocus=X<CR>*” where X could be as positive as negative values.

Returns: “*Focus=Y*”, where Y is current value of Focus motor or “*errorFocus*” when focus value is not reachable;

Action: Move Aperture motor on X step (one-quarter-stop f-number)

Command: “*moveAper=X<CR>*” where X could be as positive as negative values. **Returns:** “*Iris=Y*”, where Y is current value of Aperture or “*errorAperLimits*” when aperture value is not reachable

Action: Set desired value for Focus motor

Command: “*setFocus=X<CR>*” where X is positive value

Returns (in two UDP packets): “*Focus =XXX*”, where XXX is the current value (for non-USM lens could take more time), and “*focusDone*” message, “*errorFocus*” - indicate focus positioning problem

Action: Timeout for *fRange/setFocus/setFocus2* for old and slow lens without unreliable feedback from encoder. Timeout will be stored in non-volatile device memory.

Command: "*setFocusT=X<CR>*" where X is a positive value, presents timeout in ms. The default value is 300 ms, available range from 0 to 65000.

Returns: "OK"

Action: an alternative way to set lens focus: lens will move the focus motor to minimum position, reset the encoder counter, and then position the focus. Helps to suppress the accumulative lens drift effect for some lenses.

Command: "*setFocus2=X<CR>*" where X is a positive value

Returns(in two UDP packets): "*Focus =XXX*", where XXX is the current value (for non-USM lens could take more time), and "*focusDone*", "*errorFocus*" - indicate focus positioning problem

Action: Set desired value for aperture (in f-number)

Command: "*setAper=X<CR>*" where X is in quarter-stop f-number scale, please use pre-calculated f-stop numbers: 1.0, 1.1, 1.3, 1.4, 1.6, 1.8, 2.0, 2.2, 2.5, 2.8, 3.2, 3.5, 4.0, 4.5, 5.0, 5.6, 6.3, 7.1, 8.0, 9.0, 10, 11, 13, 14, 16, 18, 20, 22, 25, 29, 32, 36, 40, 45, 51, 57, 64, 72, 80, 90

Returns: "*Iris=Y*", where Y – current value of Aperture, if X is out of range will answer with "*errorAperLimits*" message

Action: detect IS (Image Stabilization) availability for connected lens

Command: "*isIS<CR>*"

Returns: “IS=Y” where Y could be **0** or **1**, “**0**” - means the lens does not have IS function, “**1**” - lens has IS function. Also will answer “**0**” if the IS switch is in OFF position.

Action: activate IS (Image Stabilization) function for connected lens

Command: “enableIS=X<CR>”

Returns: “ISactive=Y” where X is the number of seconds to keep IS active, from [0-3600]. “**0**” is disable IS, and Y is **0** or **1** is current IS status

Action: Set nickname for LC (stored in LC memory)

Command: “setName=XXXXXXX<CR>” where XXXXXXX is a 7-characters name for this controller.

Returns: “OK”

Action: Get LC nickname (stored in non-volatile memory)

Command: “getName<CR>”

Returns: “Name=XXXXXX”

Action: soft restart LC-2

Command: “Reboot<CR>”

Returns: “OK”

Action: explore focus limits

Command: “refRange<CR>”

Returns: “fRange=minY,maxY” where minY is minimum focus value, maxY is maximum focus value

Action: get lens name

Command: "getLens<CR>"

Returns: "Lens=XXXXXXX", where XXXXXXXX is the lens name stored in internal lens memory, the lens might not support this command.

Action: get distance range from internal lens encoder in meters

Command: "getDist<CR>"

Returns: "Dist=XX,YY", where XX and YY means that current lens focus distance is in the range of XX-YY meters

Action: set focus ring in minimum position

Command: "goMIN<CR>"

Returns: "MIN_OK"

Action: set focus ring in maximum (infinity) position

Command: "goMAX<CR>"

Returns: "MAX_OK"

Action: open aperture to the maximum for that lens and current zoom position

Command: "openAper<CR>"

Returns: "OK"

Action: reset internal encoder counter

Command: "setZero<CR>"

Returns: "OK"

Action: enable a workaround for Canon EFS 18-135mm lens

Command: “*setProto=X<CR>*”, where X is 1 or 0, will be stored in non-volatile device memory.

Returns: “OK”

Action: get current status for workaround for Canon EFS 18-135mm lens

Command: “*getProto<CR>*”,

Returns: “1” or “0”

Action: set desired lens value for zoom in mm

Command: “*setZoom=X<CR>*”, where X is lean focal lengths value in mm.

Returns(in two UDP packets): “Zoom=XXX”, where XXX is the current zoom value and “ZoomDone” message

Action: get current lens zoom value in mm

Command: “*getZoom<CR>*”,

Returns: “Zoom=XXX” where XXX is the zoom length value in mm.

Action: move zoom motor

Command: “*moveZoom=X<CR>*”, where X internal zoom steps (for precise move zoom ring), could be negative or positive.

Returns: “Zoom=XXX”, where XXX is the current zoom value

2. Examples

Python

```
import socket
import time

UDP_IP = "192.168.2.251"
UDP_PORT = 1339

counter = 1
print ("ISSI :: Canon Python Script Example 2016-8-25\n")

# LC command sequence
MESSAGES = ["ver", "ping", "moveAper=3","moveAper=-3", "Restart"]

sock = socket.socket(socket.AF_INET, socket.SOCK_DGRAM) # UDP
sock.setblocking(0)
sock.settimeout(2.5)

# main loop
while True:
    print("#%d cycle\n" %counter)
    for i in range(len(MESSAGES)): # LC2 command loop
        print ("Sent:", MESSAGES[i])
        sock.sendto(MESSAGES[i].encode(), (UDP_IP, UDP_PORT))
        try:
            while True:
                data, addr = sock.recvfrom(1024)
                if not data: break
                print ("Received:", data)
        except socket.error:
            print("")
        time.sleep(1.0)
        counter = counter +1
```

Script reply for Canon EF 35mm f/1.4L USM lens

ISSI :: Canon Python Script Example 2016-8-25

#1 cycle

Sent: ver

Received: b'ISSI LC-2=1.8.0 (s/n:19462)'

Sent: ping

Received: b'zRange=35,35'

Received: b'fRange=0,2505'

Received: b'aRange=1.4,22.0'

Received: b'Current=35,2505,1.4'

Received: b'AF=1'

Received: b'IS=0'

Received: b'ISactive=0'

Sent: moveAper=3

Received: b'Iris=2.0'

Sent: moveAper=-3

Received: b'Iris=1.4'

Sent: Restart

Received: b'OK'

[...]
